

## SPENCER 100

### Lengthwise foldable stretcher



Foldable stretcher with ground supports for patient recovery and transfer to a main stretcher or other transport device.

#### Specific features

- Extremely lightweight and durable structure
- Very low storage space required
- Lengthwise foldable
- Easy to clean sheet
- The handles offer a firm grip

#### Standard equipment

ST00598A STX598 RESTRAINT BELTS

#### Reference standards

NA

#### Technical data

Length	2010 ± 10 mm
Length folded	1000 ± 10 mm
Width	480 ± 5 mm
Ground clearance	140 ± 5 mm
Thickness folded	170 ± 5 mm
Materials	Al, Steel, PVC
Loading capacity	170 kg
Weight	5,5 ± 0,3 kg

Class I Med.Device according to Dir. 93/42/CEE

SPENCER 100 – FOLDABLE  
EMERGENCY STRETCHER YELLOW  
ST30099A



CND Classification	V9099
Registration number	172333

Rev.0 (02/10/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[www.spencer.it](http://www.spencer.it)