

SPENCER 100

Lengthwise foldable stretcher



Spencer 100 field stretchers are devices for lifting and carrying people which are also suitable for parking in field rescue conditions. If they are not equipped with supports, they can be parked using the Porta Stock stand.

Specific features

- Extremely lightweight and durable structure
- Very low storage space required
- Lengthwise foldable
- Easy to clean sheet
- The handles offer a firm grip

Technical data

Length when open (mm)	2010 ± 10 mm
Length when folded (mm)	1000 ± 10 mm
Width when open (mm)	480 ± 5 mm
Height when open (mm)	140 ± 5 mm
Height when folded (mm)	170 ± 5 mm
Weight (kg)	5.5 ± 0.3 kg
Capacity (kg)	170
Gripping points (no.)	4
Support feet (no.)	4
Materials	Al, Steel, PVC

Standard equipment

STX 598 Restraint belts 2 pcs

Class I MD compliant with UE Reg. 2017/745

SPENCER 100 FOLDABLE (LENGTH)
STRETCHER GREEN

ST30108A

CND Classification



Registration number

V9099

172336

Rev.0 (03/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it