

# SPENCER 110

## Widthwise foldable stretcher



Foldable stretcher with ground supports for patient recovery and transfer to a main stretcher or other transport device.

### Specific features

Extremely lightweight and durable structure  
Very low storage space required  
Widthwise foldable  
Easy to clean sheet  
The handles offer a firm grip

### Technical data

Length	2010 ± 10 mm
Width	550 ± 5 mm
Width folded	110 ± 20 mm
Ground clearance	140 ± 5 mm
Materials	Al, Steel,PVC
Loading capacity	170 kg
Weight	6 ± 0,3 kg

### Standard equipment

ST00598A STX598 RESTRAINT BELTS

### Reference standard

NA

Class I Med.Device according to Dir. 93/42/CEE

SPENCER 110 FOLDABLE POLE  
STRETCHER

ST30112A

CND Classification

Registration number



V9099

172343

Rev.0 (10/10/2019)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[spencer.it](http://spencer.it)