


SPENCER 110

Widthwise foldable stretcher



Spencer 110 field stretchers are devices for lifting and carrying people which are also suitable for parking in field rescue conditions. If they are not equipped with supports, they can be parked using the Porta Stock stand.

Specific features

- Extremely lightweight and durable structure
- Very low storage space required
- Widthwise foldable
- Easy to clean sheet
- The handles offer a firm grip

Technical data

Length when open (mm)	2010 ± 10 mm
Length when folded (mm)	2010 ± 10 mm
Width when open (mm)	550 ± 5 mm
Width when folded (mm)	110
Height when open (mm)	140
Height when folded (mm)	140 ± 5 mm
Weight (kg)	6 ± 0.3 kg
Capacity (kg)	170
Gripping points (no.)	4
Support feet (no.)	4
Materials	Al, Steel, PVC

Standard equipment

STX 598 Restraint belts 2 pcs

Class I MD compliant with UE Reg. 2017/745

SPENCER 110 FOLDABLE (WIDTH)
STRETCHER GREEN
ST30113A



CND Classification	V9099
Registration number	172346

Rev.0 (03/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
spencer.it