


## NRT

### Recovery stretcher



NRT recovery stretchers are devices to be used when the morphological characteristics of the intervention area require particularly compact equipment and adherence to the patient.

#### Specific features

- Low storage space required
- Transport handles allow easy use
- Integrated adjustable belts for patient
- Belt for head stabilization

#### Technical data

Dimensions	1080 x 1830 $\pm$ 5 mm
Thickness	20 $\pm$ 5 mm
Material	Wood, Cotton
Number of handles	6
Loading capacity	220 kg (Using all handles simultaneously)
Weight	7 $\pm$ 0,5 kg

Class I MD compliant with UE Reg. 2017/745

NRT – WHITE RECOVERY  
STRETCHER

**ST30150A**

**CND Classification** V9099

**Identification number** 216053



Rev.0 (24/05/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[spencer.it](http://spencer.it)