

## SPENCER 200

### Foldable stretcher with wheels



Spencer 200 field stretchers are devices for lifting and carrying people which are also suitable for parking in field rescue conditions. If they are not equipped with supports, they can be parked using the Porta Stock stand.

#### Specific features

- Extremely lightweight and durable structure
- Very low storage space required
- Lengthwise foldable
- Easy to clean sheet
- Carrying handles offer a firm grip
- 2 wheels  $\varnothing$  100 mm
- 2 fixed support feet

#### Technical data

Length when open (mm)	1885 $\pm$ 5 mm
Length when folded (mm)	925 $\pm$ 5 mm
Width (mm)	490 $\pm$ 5 mm
Thickness when folded (mm)	165 $\pm$ 5 mm
Weight (kg)	5 $\pm$ 0.2 kg
Maximum load capacity (kg)	170
Wheels	2 each 100 mm
Support feet	2
Materials	Al, Steel, PVC

#### Standard equipment

STX 598 Restraint belts 2 pcs

Class I MD compliant with UE Reg. 2017/745

SPENCER 200 EMERGENCY  
FOLDING STRETCHER ORANGE  
ST30200A



CND Classification	V9099
Registration number	168576

Rev.0 (03/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[spencer.it](http://spencer.it)