


SPENCER 250

Lengthwise compact foldable stretcher



Spencer 250 field stretchers are devices for lifting and carrying people which are also suitable for parking in field rescue conditions. If they are not equipped with supports, they can be parked using the Porta Stock stand.

Specific features

- Extremely lightweight and durable structure
- Very low storage space required
- Lengthwise foldable
- Easy to clean sheet
- The high number of handles make transport easier

Technical data

Length when open (mm)	1890 ± 10 mm
Length when folded (mm)	950 ± 10 mm
Width (mm)	480 ± 5 mm
Thickness (mm)	30 mm
Thickness when folded (mm)	60 ± 5 mm
Weight (kg)	5.5 ± 0.3 kg
Maximum load capacity (kg)	170
Materials	Al, PVC

Standard equipment

STX 598 Restraint belts 2 pcs

Class I MD compliant with UE Reg. 2017/745

EASY-FOLD EMERGENCY
STRETCHER SPENCER 250
ST30250A



CND Classification	V9099
Registration number	168568

Rev.0 (03/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it