

SPENCER 270

Emergency stackable stretcher



Spencer 270 field stretchers are devices for lifting and carrying people which are also suitable for parking in field rescue conditions. If they are not equipped with supports, they can be parked using the Porta Stock stand.

Specific features

- Extremely lightweight and durable structure
- Low storage space required thanks to the possibility to stack many stretchers of the same model
- Easy to clean PVC sheet
- 8 handles and the dimensions of the tubular frame allow a firm grip during transfer procedures
- 4 support feet confer high stability

Technical data

Length when open (mm)	1885 ± 10 mm
Width (mm)	480 ± 5 mm
Height (mm)	175 ± 10 mm
Weight (kg)	5 ± 0.5 kg
Capacity (kg)	170
Gripping points (no.)	8
Support feet (no.)	4
Approximate overall dimensions 10 stacked stretchers	610 ± 30 mm
Materials	Al, Steel, PVC

Standard equipment

- STX 598 Restraint belts 2 pcs

Class I MD compliant with UE Reg. 2017/745

SPENCER EMERGENCY STRETCHER
MOD. 270

ST30270A

CND Classification V9099

Registration number 94048



Rev.0 (03/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it