

SPENCER 280

Emergency stackable stretcher



Spencer 280 field stretchers are devices for lifting and carrying people which are also suitable for parking in field rescue conditions. If they are not equipped with supports, they can be parked using the Porta Stock stand.

Specific features

- Extremely lightweight and durable structure
- Low storage space required thanks to the possibility to stack many stretchers of the same model
- Easy to clean PVC sheet
- The dimensions of the tubular frame allow a firm grip during transfer procedures
- The support feet confer high stability

Technical data

Length when open (mm)	1885 ± 10 mm
Width (mm)	480 ± 5 mm
Height (mm)	175 ± 10 mm
Weight (kg)	6 ± 0.5 kg
Capacity (kg)	170
Gripping points (no.)	8
Support feet (no.)	4
Materials	Al, Steel, PVC

Standard equipment

- STX 598 Restraint belts 2 pcs

Class I MD compliant with UE Reg. 2017/745

STACKABLE STRETCHER SPENCER
280

ST30280A

CND Classification V9099

Registration number 94039



Rev.0 (03/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it