

SPENCER 281

Emergency stackable stretcher with head support



Spencer 281 field stretchers are devices for lifting and carrying people which are also suitable for parking in field rescue conditions. If they are not equipped with supports, they can be parked using the Porta Stock stand.

Specific features

- Extremely lightweight and durable structure
- Very low storage space required thanks to the possibility to stack other stretcher of the same model
- Backrest adjustable in 4 position other than the horizontal one
- PVC sheet Easy to clean
- The dimensions of the tubular frame allow a firm grip during transfer procedures
- The support feet confer high stability

Technical data

Length when open (mm)	1885 ± 10 mm
Width (mm)	480 ± 5 mm
Height (mm)	175 ± 10 mm
Weight (kg)	6 ± 0.5 kg
Capacity (kg)	170
Gripping points (no.)	6
Support feet (no.)	4
Materials	Al, Steel, PVC

Standard equipment

- STX 598 Restraint belts 2 pcs

Class I MD compliant with UE Reg. 2017/745

STOCKABLE EMERGENCY
STRETCHER WITH BACKREST
ORANGE



ST30281A

CND Classification	V9099
Registration number	94200

Rev.0 (03/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it