


SPENCER 310

Foldable stretcher with backrest



Spencer 310 field stretchers are devices for lifting and carrying people which are also suitable for parking in field rescue conditions. If they are not equipped with supports, they can be parked using the Porta Stock stand.

Specific features

- Extremely lightweight and durable structure
- Very low storage space required
- Lengthwise foldable
- Easy to clean sheet
- Adjustable backrest increase patient comfort

Technical data

Length when open (mm)	1870 ± 10 mm
Length when folded (mm)	945 ± 10 mm
Width (mm)	500 ± 5 mm
Thickness when folded (mm)	160 ± 5 mm
Height from the ground	160 ± 5 mm
Weight (kg)	6.5 ± 0.3 kg
Maximum load capacity (kg)	170
Wheels	2 each 100 mm
Support feet	2
Materials	Al, PVC

Standard equipment

STX 598 Restraint belts 2 pcs

Class I MD compliant with UE Reg. 2017/745

SPENCER 310 EMERGENCY
FOLDING STRETCHER WITH
BACKREST ORANGE



ST30310

CND Classification V9099

Registration number 168589

Rev.0 (03/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it