


SPENCER 310

Foldable stretcher with backrest



Stretcher for recovery and transfer of patients to a main stretcher or other carrying devices.

Specific characteristics

- Extremely lightweight and durable structure
- Very low storage space required
- Lengthwise foldable
- Easy to clean sheet
- Adjustable backrest increase patient comfort

Technical data

Length	1870 ± 10 mm
Length folded	945 ± 10 mm
Width	500 ± 5 mm
Height form ground	160 ± 5 mm
Thickness folded	160 ± 5 mm
Wheels	Ø100 mm
Backrest positions	4
Materials	Al,PVC
Weight	6,5 ± 0,3 kg

Class I Med.Device according to Dir. 93/42/CEE

MOD.E020 SPENCER 310 FOLD.

STRETCHER W/ BACKREST

ST30312

CND Classification



V9099

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
www.spencer.it