

Lower fastener

For chairs 420/425/450/455



The fastener is intended for installation inside a medical vehicle. The device is compatible with Spencer evacuation chairs.

Specific Characteristics

Compact and robust
Foot operated simple and intuitive mechanism

Technical data

Base dimensions	401x72x80 ± 5 mm
Material	Steel
Weight	1 ± 0,1 kg

Class I Med.Device according to Dir. 93/42/CEE

MOUNTING BRACKET FOR CHAIR
SPENCER 450-55/420-25

ST30451B

Classification CND V9099

Registration number 1502133



Rev.0 (16/11/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it