

Spencer 456 Avio

Transport chair with four wheels



Spencer 456 Avio is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle. The reduced width allows it to be used in confined spaces.

Specific Characteristics

- Easy to clean PVC sheets
- 4 rear foldable handles to reduce encumbrance and to facilitate the lifting in every situation
- Telescopic front handles to facilitate the handling
- Foldable to save space
- Strong and light frame
- Foldable footrest to increase stability and comfort of the patient
- Front wheels with brake for rotation and axis

Standard equipment

- ST30463A Pair of belts for chair 451/456

Accessories

- ST40428B Stabilizer stuffed seat/backrest profile series

Technical data

Length opened w/closed handles	610 ± 10 mm
Length opened w/opened handles	1140 ± 10 mm
Width	450 ± 10 mm
Height opened	920 ± 10 mm
Thickness folded	235 ± 10 mm
Rear wheels	Ø125 mm
Front wheels	Ø100 mm
Materials	Acciaio/Al/PU/PVC
Weight	11,6 kg ± 0,3
Loading capacity	180 kg

Class I Med.Device according to Dir. 93/42/CEE

SPENCER CHAIR 456 AVIO W/ 4

WHEELS

ST30456A

CND Classification

Registration number



V9099

93484

Rev.3 (28/01/202)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
spencer.it

©Copyright - This document is property of SPENCER ITALIA S.r.l and can't be reproduced, used and disclosed without the written approval of SPENCER ITALIA S.r.l.