

## Spencer 456 Avio

### Transport chair with four wheels



Spencer 456 Avio is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle. The reduced width allows it to be used in confined spaces.

#### Specific Characteristics

Easy to clean PVC sheets
4 rear foldable handles to reduce encumbrance and to facilitate the lifting in every situation
Telescopic front handles to facilitate the handling
Foldable to save space
Strong and light frame
Foldable footrest to increase stability and comfort of the patient
Front wheels with brake for rotation and axis

#### Standard equipment

ST30463A	Pair of belts for chair 451/456
----------	---------------------------------

#### Accessories

ST36458B	F 400-A fixation system 20 g for chairs
ST40428B	Stabilizer stuffed seat/backrest profile series

#### Technical data

Length opened w/closed handles	610 ± 10 mm
Length opened w/opened handles	1140 ± 10 mm
Width	450 ± 10 mm
Height opened	920 ± 10 mm
Thickness folded	235 ± 10 mm
Rear wheels	ø125 mm
Front wheels	ø100 mm
Materials	Acciaio/Al/PU/PVC
Weight	11,6 kg ± 0,3
Loading capacity	180 kg

#### Class I Med.Device according to Dir. 93/42/CEE

SPENCER CHAIR 456 AVIO W/ 4 WHEELS

ST30456A

CND Classification V9099

Registration number 93484



Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[www.spencer.it](http://www.spencer.it)