

## Spencer 458/B

### 10G Transport chair



Spencer 458B is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

#### Specific Characteristics

Easy to clean PVC sheets
4 rear foldable handles to reduce encumbrance and to facilitate the lifting in every situation
Telescopic front handles to facilitate the handling
Foldable to save space
Strong and light frame
Foldable footrest to increase stability and comfort of the patient
Front wheels with brake for rotation and axis
Foot operated centralized brake for rear wheels

#### Technical data

Length opened w/closed handles	690 ± 10 mm
Length opened w/opened handles	960 ± 10 mm
Width	510 ± 10 mm
Height opened	910 ± 10 mm
Thickness folded	260 ± 10 mm
Rear wheels	Ø150 mm
Front wheels	Ø100 mm
Materials	Steel/Al/PU/PVC
Weight	12 kg ± 0,3
Loading capacity	180 kg

#### Standard equipment

ST30443A	Pair of belts for chair 450/455
----------	---------------------------------

#### Accessories

ST36456B	C-Max 10g Fixing System <small>Mandatory for 10g versions</small>
ST40428B	Stabilizer stuffed seat/backrest profile series

Class I Med.Device according to Dir. Regulation 2017/745/UE

TRANSPORT CHAIR MOD.458  
CERTIFIED 10G

**ST36453B**

**CND Classification**



V9099