

## Spencer 458/B

### 10G Transport chair



Spencer 458B is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

#### Specific Characteristics

- Easy to clean PVC sheets
- 4 rear foldable handles to reduce encumbrance and to facilitate the lifting in every situation
- Telescopic front handles to facilitate the handling
- Foldable to save space
- Strong and light frame
- Foldable footrest to increase stability and comfort of the patient
- Front wheels with brake for rotation and axis
- Foot operated centralized brake for rear wheels
- Padded armrest to improve patient containment and comfort

#### Standard equipment

ST30443A	Pair of belts for chair 450/455
----------	---------------------------------

#### Accessories

ST36456B	C-Max 10g Fixing System <small>Mandatory for 10g versions</small>
ST40428B	Stabilizer stuffed seat/backrest profile series

#### Technical data

Length opened w/closed handles	690 ± 10 mm
Length opened w/opened handles	960 ± 10 mm
Width	510 ± 10 mm
Height opened	910 ± 10 mm
Thickness folded	260 ± 10 mm
Rear wheels	Ø150 mm
Front wheels	Ø100 mm
Materials	Steel/Al/PU/PVC
Weight	12 kg ± 0,3
Loading capacity	180 kg

#### Class I Med.Device according to Dir. 93/42/CEE

TRANSPORT CHAIR MOD.458/B

CERTIFIED 10G

ST36454B

CND Classification

Registration number



V9099

93458