

## Spencer 455

### 20G Transport chair with four wheels



Spencer 455 is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

#### Specific Characteristics

Easy to clean PVC sheets
4 rear foldable handles to reduce encumbrance and to facilitate the lifting in every situation
Telescopic front handles to facilitate the handling
Foldable to save space
Strong and light frame
Foldable footrest to increase stability and comfort of the patient
Front wheels with brake for rotation and axis

#### Standard equipment

ST30443A	Pair of belts for chair 450/455
----------	---------------------------------

#### Accessories

ST36451B	F 400 fixation system 20 g for chairs
ST40428B	Stabilizer stuffed seat/backrest profile series

#### Technical data

Length opened w/closed handles	610 ± 10 mm
Length opened w/opened handles	1140 ± 10 mm
Width	510 ± 10 mm
Height opened	920 ± 10 mm
Thickness folded	250 ± 10 mm
Rear wheels	ø125 mm
Front wheels	ø100 mm
Materials	Steel/Al/PU/PVC
Weight	12 kg ± 0,3
Loading capacity	180 kg

Class I Med.Device according to Dir. 93/42/CEE

STAIR CHAIR SPENCER 455 WITH 4  
WHEELS  
ST36455A



CND Classification	V9099
Registration number	413100