

Roller AV

Self-Loading trolley with adjustable heights



Self-loading trolley intended for transport and support of Spencer stretchers series Roller.

Specific features

Extremely versatile thanks to the adjustable heights system
 Leg's closure commands can be easily reached
 Swiveling rear wheels
 Brake on front wheels
 The button to unlock the stretcher from the trolley can be easily reached
 The size and the material of the wheels minimizes friction on smooth surfaces

Technical data

Length	1900 mm
width	560 ± 10mm
Rised trolley height ¹	820 ± 10mm
Folded trolley height	230 ± 10mm
Loading trolley wheel height	69 cm
Possible heights	Raised Intermediate 1 Intermediate 2 Folded
Materials	Al, Steel, PE
Loading capacity	170 kg
Weight	38 ± 1 kg

¹Dimensions referred to the stretcher support area

Optionals

ST42100A	Rear fastener FP
ST42200A	Front fastener FA

Class I Med.Device according to Dir. 93/42/CEE

ROLLER AV TROLLEY IN 3 LEVELS
 LOADING HEIGHT 69

ST40019A

CND Classification



Registration number

V0899

159685

Rev.0 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it