

# ROLLER S

## Stretcher with telescopic handles



Optional mattress

Stretcher for recovery of patients to be used in cases where the use of monoblock stretchers is not possible or difficult. The device can be used with Spencer Roller AC/AV self-loading trolleys.

### Specific features

- Robust and lightweight aluminium frame
- The wheels allow easy movements on the recovery area
- Adjustable backrest
- The use with the dedicated trolley increases its versatility as well as safety during transport
- The attachment on the trolley doesn't require additional devices
- Foldable side bars *Safe Bar* increase patient safety

### Technical data

Length w/handles closed	1930 ± 10 mm
Length w/handles opened	2340 ± 10 mm
Width	570 ± 10 mm
Wheels	∅ 100 mm
Side bars dimensions	680x200 mm
Materials	Al, PE, Steel
Weight	22 ± 1 kg
Loading capacity	170 kg

### Standard equipment

ST00526A	STX526 2 PCS METAL ORANGE BELTS W/METAL BUCKLE
----------	--

### Optionals

ST70000A	QMX/777 PVC anatomic mattress
----------	-------------------------------

Class I Med.Device according to Dir. 93/42/CEE

ROLLER S ADJ.BRACKETS  
TELESCOPIC HANDLES



ST40300A

CND Classification	V0899
Registration number	159689

Rev.0 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)