

## Spencer 407/B

Transport chair with 4 wheels



Spencer 407/B is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

### Specific features

Easy to clean PVC sheets
Armrest and rear handles made of a single component, allowing weight reduction
Robust and lightweight structure
Footrest and armrests to improve the stability and patient comfort
Large rear wheels to reduce vibration during patient transfer
Swiveling front wheels with brake for rotation and axis
Quick and intuitive closure system
Safety system against accidental closure

### Reference standard

NA

### Standard equipment

ST10408A | 3 belts for chair 406/407

### Optional accessories

ST42705E | 10g wall fastener for chair 407

### Technical data

Length with handles closed	810 ± 10 mm
Length with handles opened	880 ± 10 mm
Width	475 ± 10 mm
Height opened	970 ± 10 mm
Height folded	800 ± 10 mm
Thickness folded	250 ± 10 mm
Rear wheels	∅ 200 mm
Front wheels	∅ 100 mm
Materials	Steel/Al/PVC
Weight	12,5 kg ± 0,3
Loading capacity	150 kg

Class I Med.Device according to Dir. 93/42/CEE

SPENCER 407B CHROME  
W/ARMRESTS EVACUATION CHAIR  
4 WHEELS



**ST40407A**

**CND Classification** V9099

**Registration number** 93563

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)