

Spencer 455/B

Transport chair with four wheels and armrests



Spencer 455/B is a transport chair for the evacuation and transport of patient from the rescue place to the ambulance or to another device intended for transport on the vehicle.

Specific Characteristics

Easy to clean PVC sheets
4 rear foldable handles to reduce encumbrance and to facilitate the lifting in every situation
Telescopic front handles to facilitate the handling
Foldable to save space
Strong and light frame
Foldable footrest to increase stability and comfort of the patient
Front wheels with brake for rotation and axis

Standard equipment

ST30443A	Pair of belts for chair 450/455
----------	---------------------------------

Accessories

ST30451A	MOUNTING BRACKETS FOR CHAIR SPENCER 450-55/420-25
ST36452E	UPPER FIXATION SYSTEM FOR TRANSPORT CHAIRS
ST40428B	Stabilizer stuffed seat/backrest profile series

Reference standard

NA

Technical data

Length opened w/closed handles	610 ± 10 mm
Length opened w/opened handles	1140 ± 10 mm
Width	510 ± 10 mm
Height opened	920 ± 10 mm
Thickness folded	250 ± 10 mm
Rear wheels	ø125 mm
Front wheels	ø100 mm
Materials	Steel/Al/PU/PVC
Weight	12,7 kg ± 0,3
Loading capacity	180 kg

Class I Med.Device according to Dir. 93/42/CEE

SPENCER 455/B WITH 4 WHEELS
AND ARMRESTS



ST40455A

CND Classification V9099

Registration number 93475