

FP

Rear fastener



System to be used to fix the stretchers on the floor of the medical vehicle. The device is to be used together with a front fastener. This version is specific for stretchers series 506.

Specific features

- Dimensions reduced to the minimum
- Easy locking and unlocking
- Requires a small surface for its installation

Technical data

Length	131 mm
Width	115 mm
Height with high lock	161 mm
Height with low lock	114 mm
Materials	Cast iron
Weight	1,8 ± 0,1 kg

Class I Med.Device according to Dir. 93/42/CEE

FA ANTERIOR FASTNER FOR SELF
LOADING STRETCHERS

ST42100A



CND Classification V9099

Registration number 91964

Rev.1 (23/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it