


FA - Front fastener for stretcher



Fastener for stretchers, created to provide an anchorage inside medical vehicles. The device should be used together with a rear fastener.

Specific features

- The compact structure avoids unnecessary bulk when the stretcher is not loaded on the ambulance
- Allows easy installation
- The weight is minimized

Technical data

Dimensions	455x430x67 ± 5mm
Materials	Stainless Steel

Class I Med.Device according to Dir. 93/42/CEE

FRONT FASTENER FOR
STRETCHERS
ST42200A

CND Classification



V9099

Rev.1 (23/05/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it