


## R MAX

### 10G fixing system for stretcher



Fixing system for 10g certified stretchers, created to ensure their safe anchorage inside medical vehicles.

#### Specific Characteristics

Reduced encumbrance  
Coupling mechanism simple and intuitive  
3 fixing point for the stretcher for a better distribution of forces

#### Compliant with

EN 1789

#### Technical data

Dimensions	1775x445x75 ± 15mm
Weight	9,6 ± 0,2 kg

#### Class I Med.Device according to Dir. 93/42/CEE

R-MAX B 10G HOOK FOR  
STRETCHERS  
ST42702C



CND Classification	V9099
Registration number	92046

Rev.3 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)