

10G Wall Fastener



10g wall fastener to be installed inside an ambulance, suitable for transport chairs Spencer 406 and 407.

Specific features

Compact and robust
 Easy to install, only 3 installation holes are required
 Extremely easy to use
 Does not require additional components to mount the chair

Technical data

Wall mounted dimensions	645x305x150 mm
Materials	Al, Nylon, Steel
Weight	1,65 ± 0,1 kg

Reference standard

EN 1789

Class I Med.Device according to Dir. 93/42/CEE

10G WALL FASTENER FOR CHAIR
 406/407

ST42705E

CND Classification



V9099

Rev.2 (09/08/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it