

SENSOR LOCK

10G Fastener for stretcher



Fastener for 10g certified stretchers, created to ensure their safe anchorage inside medical vehicles. The device integrates an innovative electromechanical unlocking system that, by interfacing with the Sharp Sensor Lock system installed on the stretcher, allows its release without having to act on the fixing system.

Specific features

The linearity of the structure avoids unnecessary bulk when the stretcher is not loaded on the ambulance

Semiautomatic release of the stretcher

The device acts simultaneously both on frontal and rear part of the stretcher and, thanks to the integral structure, allows a simple installation without alignment difficulties between the frontal and rear part of the fastener

Equipped with advanced maintenance control systems

Increased safety for operators and patients thanks to information signals to show if the stretcher is fastened/not fastened/not properly fastened

The weight is minimized thanks to researches of advanced materials and in-depth studies about constructive geometries

Technical data

Dimensions	1790x190x80 ± 5mm
Metal profile thickness	30 mm
Materials	Al , Nylon, PE, Steel

Compliant with

EN 1789

Class I Med.Device according to Dir. 93/42/CEE

10G FASTENER SENSOR LOCK

ST42706



CND Classification

V9099

Rev.0 (17/04/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it