


## FA HUX

### Anterior Fastener for stretcher



Anterior fastener for stretchers, created to provide an anchorage inside medical vehicles.

#### Specific features

The linearity of the structure avoids unnecessary bulk when the stretcher is not loaded on the ambulance

The release of the stretcher happens easily, simply operating on the small release mechanism on the rear fastener

The weight is minimized thanks to researches of advanced materials and in-depth studies about constructive geometries

#### Technical data

Dimensions	1625x190x80 ± 20mm
Metal profile thickness	30 mm
Materials	Al , Nylon, PE, Steel
Weight	3,8 ± 0,2 kg

Class I Med.Device according to Dir. 93/42/CEE

FA-HUX FRONT FIXATION FOR  
STRETCHER

ST42710A

CND Classification



V9099

Rev.0 (25/09/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)