

## QMX 710

### Articulated mattress for stretcher



Mattress intended for use with stretchers, in order to increase patient comfort by absorbing the stresses to which it may be subjected during the transfer stages.

#### Specific features

Smooth and easy to clean surface

The four side zips allow easy access to the padding for a possible replacement

The 4-segment body makes it particularly suitable for use with stretchers equipped with Trendelenburg / Fowler-type floor plates

Strap insert on the lower side to increase stability when placed on the stretcher

#### Technical data

Dimensions	1820x390x50 ± 10 mm
Materials	PVC, PU esp.
Strap area dimensions	250x100 ± 10 mm
Weight	2,1 ± 0,1 kg

#### Accessories

DR60851K | DRX 850 BOX 100 PCS - TNT SHEET

Class I MD. According to Dir. 93/42/CEE

QMX 710 ARTICULATED MATTRESS  
FOR STRETCHERS ORANGE

ST43010A

CND Classification V9099

Registration number 93114



Rev.2 (18/07/2017)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>

Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)