

SPENCER 506

Stretcher with telescopic handles



Ella series field stretchers Spencer 506 are devices for lifting and carrying people for use in field rescue situations. The foldable, lockable legs allow for transferring either by means of wheels on the legs or by lifting when the legs are locked in the closed position. Parking is permitted with the stretcher alone or with a trolley with closed legs.

Specific features

- Sturdy and lightweight aluminium frame
- The wheels allow easy movements in the rescue area
- Adjustable backrest
- Trendelenburg antishock position
- High versatility and safety during transport if used with the dedicated trolley
- The mounting on the trolley doesn't require additional equipment
- The side bars *Safe bar* increase patient safety

Technical data

Length with handles closed	1930 ± 10 mm
Length with handles open	2340 ± 10 mm
Width	570 ± 10 mm
Wheels	∅ 100 mm
Dimensions sidebars	680 x 200 mm
Materials	Al, PE, steel
Weight	22 ± 1 kg
Capacity	170 kg

Standard equipment

- STX526 Reflective restraint belts 2 pieces with metal buckle

Accessories

ST70000A	QMX 777 Black PVC anatomic mattress
ST10602C	510 Self trolley silver for 506 T

Class I MD compliant with UE Reg. 2017/745

506S SILVER STRETCHER
W/TELES. HANDLES FOR 510
TRO



ST60505A

CND Classification	V0899
Registration number	165270

Rev.0 (03/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it