

STX 702

Metal buckle belt



STX 702 restraint belts are aids used to obtain a degree of restraint that allows patients to be transported safely when used in accordance with the product in use.


Specific features

- Metal buckle to provide more safety
- Reflective insert inside the fabric to increase visibility
- Suitable for self-loading and basket stretchers

Technical data ¹

	Min (mm)	Max (mm)
Length L	600	1600
Tape width	50 mm	
Materials	Nylon, chrome-plated steel, plastic	
Weight of the pair	460 g	
Belt type	2 pieces	
Hook	Metal	

¹ The indicated sizes are in reference to the maximum dimensions. The useful length can be different, due to the device ore to the knot used to attach the belt to the frame of the stretcher.



Class I MD compliant with UE Reg. 2017/745

STX 702 2pcs BLACK REFLEX BELTS
w/METAL BUCKLE

ST70002A

CND Classification

Registration number



V9099

96045

Rev.0 (08/01/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it