

# QMX777

## Anatomic mattress



Mattress to be used with stretchers in order to increase patient comfort by absorbing vibrations and stresses that may occur during transfer procedures.

### Specific features

- Smooth and easy to clean surface
- Tape with strap elements, used to fix the mattress to the stretcher board
- The shape, the internal divisions and the filling material, allow excellent adaptation to the patient's body
- Compatible with boards type S and T

### Accessories

- DR60851K DRX 850 BOX 100 PCS - TNT SHEET

### Technical data

Dimensionis	1770x430 ± 10 mm
Central part thickness	50 mm
Maximum thickness	90 mm
Materials	PU exp, PVC,
Weight	2,3 ± 0,1 kg

### Class I MD. According to Dir. 93/42/CEE

QMX 777 ORANGE ANATOMIC MATTRESS, WATERTIGHT, PVC  
ST70004A



CND Classification	V9099
Registration number	93128

Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)