

QMX776

One segment mattress



Mattress to be used with stretchers in order to increase patient comfort by absorbing vibrations and stresses that may occur during transfer procedures.

Specific features

- Smooth surface easy to clean
- Strap tapes to fix the mattress to the stretcher
- The filling material allows an excellent adaptation to the body of the patient
- Compatible with patient floors type S and T

Technical data

Dimensions	1800X430X70 ± 10 mm
Materials	PU exp, PVC
Weight	2,1 ± 0,1 kg

Accessories

DR60851K	DRX 850 – DISPOSABLE GREEN TNT SHEET 100PCS BOX
----------	--

Class I MD. According to Dir. 93/42/CEE

QMX 776 BLACK ANATOMIC
MATTRESS, PVC, 1 SECTION
ST70018A



CND Classification	V9099
Registration number	93138

Rev.2 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it