

## PORTA STOCK

Support frame for emergency stretchers



Porta Stock is a support for field stretchers that allows the stretcher to be kept raised from the ground, facilitating small movements and access to the patient by the doctor.

Available in versions with and without wheels, it is compatible with the Stock series stretchers and is also suitable for Spencer stretchers with a maximum width of 575 mm and a minimum length of 1250 mm (not including handles).

Unsupervised parking is not permitted.

### Specific features

Robust and lightweight frame  
 Reduced thickness when folded, making it ideal for field use  
 Designed for use with Spencer stretchers Stock Series, is also suitable for use with Spencer emergency stretchers having maximum width of 575 mm and minimum length (handles excluded) of 1250 mm

### Technical data

Length when open	1325 ± 10 mm
Width	640 ± 10 mm
Height when open	750 ± 10 mm
Overall dimensions when folded	870 x 640 x 130 ± 10 mm
Weight	10 ± 0,5 kg
Capacity	150 kg
Frame material	Steel

Class I MD compliant with UE Reg. 2017/745

PORTA STOCK SUPPORTING  
 FRAME FOR EMERGENCY STRETCHERS

ST70673C



CND Classification V0899

Registration number 1241247

Rev.0 (03/08/2021)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[www.spencer.it](http://www.spencer.it)