

# FLEX MASK

## Mouth to mouth mask



Disposable face mask with pre inflated pillow, to be used as a protection between the rescuer and the patient during CPR procedures.

### Specific features

Transparent shell to improve patient monitoring  
 Inlet with cap to allow connection of an oxygen supply  
 Supplied with unidirectional valve

**Latex Free**

### Technical data

Dimensions	122x95x100 ± 2 mm
Materials	PVC
Weight	50± 5g

Class IIa Med. Dev. according to dir. 93/42/CEE

FLEX MASK – MOUTH TO MOUTH MASK W/VALVE

TA09012A

CND Classification



R0301019

Registration number

1173357

Rev.0 (28/03/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

[spencer.it](http://spencer.it)