

FLEX MASK

Mouth to mouth mask



Disposable face mask with pre inflated pillow, to be used as a protection between the rescuer and the patient during CPR procedures.

Specific features

- Transparent shell to improve patient monitoring
- Inlet with cap to allow connection of an oxygen supply
- Hard case for the best protection

Technical data

Mask dimension	122x95x100 ± 2 mm
Box dimensions	135x110x50 ± 2 mm
Materials	PVC, Nylon
Mask weight	45 ± 5 g
Overall weight	75 g

Ass. Med. Dev. According to dir. 93/42/CEE

FLEX MASK SET WITH
ORANGE PVC CASE INS. ON
THE BACK
TA09019A



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it