

FLEX MASK

Mouth to mouth mask



Disposable face mask with pre inflated pillow, to be used as a protection between the rescuer and the patient during CPR procedures. In this versions is supplied with soft case, disinfectant towel and latex gloves.

Specific features

Transparent shell to improve patient monitoring
 Inlet with cap to allow connection of an oxygen supply
 Case with zip closure and strap for belt attachment
 2 disinfectant towel and 2 latex gloves size L (not shown in the picture)

Technical data

Mask dimension	122x95x100 ± 2 mm
Case dimension	120x145x65 ± 5 mm
Materials	PVC, Nylon, Latex
Mask weight	45 ± 5 g
Kit weight	140 ± 10 g

Ass. Med. Dev. According to dir. 93/42/CEE

FLEX MASK WITH
 YELLOW/BLACK SOFT CASE
 TA09037A



Rev.1 (18/07/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

www.spencer.it