

TRI ROLL

Large soft case with trolley



Soft trolley ideal to protect and store cardiorespiratory reanimation kits, medication kits and diagnostic systems.

Specific features

- PVC fabric easy to clean
- Reflective elements on the belt placed on the front
- Handles and trolley to allow easy transport
- Zip closure and wide opening for easy access to the content
- 5 internal side pockets with zip closure
- 3 internal height adjustable belts for fixation of devices
- Trolley handle with 2 extension steps

Technical data

| | |
|--------------------------|------------------|
| Dimensions | 800x400x350 |
| Volume | 100 ± 5 l |
| Material | PVC, Nylon |
| Trolley handle extension | Approx 25 – 37cm |
| Weight | 4,2 kg |

TRI ROLL - LARGE SOFT CASE
TT02210C

Rev.1 (08/02/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

spencer.it