

Z-Straps

Belt for recovery board Z5005



Spider belt to be used with recovery board Z5005.

Specific features

Wide range of adjustments thanks to the strap closures
5 transversal bands to optimize immobilization

Technical data

Materials	Nylon, PP
Central belt length	1160 ± 20 mm

Z-STRAPS RESTRAINT SYSTEM FOR
Z5005
ZC00101A

Rev.1 (15/03/2018)

COPIA NON CONTROLLATA – successive revisioni saranno disponibili al link <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

mortuary.spencer.it