

Z501F

Mortuary trolley



Foldable mortuary trolley for support and transport of coffins.

Specific Characteristics

Aluminum sturdy and light weight structure
 Large wheels for easy movement
 2 brakes with axis and rotation locking
 Allows easy and fast opening e closure
 Extremely compact dimensions

Technical data

Length opened <small>(travel direction)</small>	1180 ± 20 mm
Width	480 ± 5 mm
Dimensions folded	800x430x130 ± 10mm
Wheels dimensions	ø200 mm
Brakes	2
Weight	12,5 ± 0,2 kg
Loading capacity	300 kg

Accessories

ZC00831A	Carry bag for trolley Z501F
ZC01850B	Extractable cross for church truck silver
ZC01854B	Extractable polished cross for church truck

Z501F SILVER CHURCH TRUCK ALU
 MINUM W/WHEELS 200M
ZC00813A

Rev.2 (15/01/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

mortuary.spencer.it