

Z2101

Recovery stretcher



Stretcher for recovery of corps.

Specific features

Resistant PVC fabric and easy to clean
Footrest and wheels on the bottom, allow to move the device by tilting and sliding it on the ground

Technical data

Length	2005 ± 10 mm
Width	485 ± 10 mm
Thickness w/footrest opened	350 ± 10 mm
Loading capacity	200 kg

Z2101 RECOVERY STRETCHER
ZS00211C

Rev.0 (25/08/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

mortuary.spencer.it