

Z2200

Foldable recovery stretcher with body bag



Foldable stretcher for recovery of corps, equipped with body bag.

Specific features

- Resistant PVC fabric and easy to clean
- Footrest and wheels on the bottom, allow to move the device by tilting and sliding it on the ground
- The cover allow covering of the corp
- The device can be folded, minimizing storage space

Technical data

Length	2005 ± 10 mm
Width	485 ± 10 mm
Thickness w/footrest opened	350 ± 10 mm
Dimensions folded	1200x485x190 ± 10 mm
Loading capacity	200 kg

Z2200 FOLDABLE RECOVERY
STRETCHER W/COVER
ZS00220C

Rev.0 (25/08/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

mortuary.spencer.it