

Z2201

Foldable recovery stretcher



Foldable stretcher for recovery of corps.

Specific features

- Resistant PVC fabric and easy to clean
- Footrest and wheels on the bottom, allow to move the device by tilting and sliding it on the ground
- The device can be folded, minimizing storage space

Technical data

| | |
|-----------------------------|----------------------|
| Length | 2005 ± 10 mm |
| Width | 485 ± 10 mm |
| Thickness w/footrest opened | 350 ± 10 mm |
| Dimensions folded | 1200x485x190 ± 10 mm |
| Loading capacity | 200 kg |

Z2201 FOLDABLE RECOVERY
STRETCHER
ZS00221C

Rev.0 (25/08/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

mortuary.spencer.it