

## Z 5000

### Pickup stretcher



The pickup stretcher Spencer Z5000 allows the operator to raise and transfer the corpse without having to move it.

#### Specific characteristics

- Extremely light and resistant structure
- The possibility to adjust the length confers high versatility
- Foldable to save space when stored
- Smooth surfaces easy to clean
- The numerous handles make the transport easier
- Belts with transport bag

#### Standard accessories

ZS00401A	ZSX401 SET 3 BELTS W/BAG
----------	--------------------------

#### Optional accessories

ST05010A	PVC storage and carry bag
ST05018C	SXV – 20G vertical fixing system
ST05019C	SXO – 20G horizontal fixing system

#### Technical data <sup>1</sup>

Shortened length	1660 mm
Intermediate length 1	1780 mm
Intermediate length 2	1900 mm
Maximum length	2050 mm
Folded length	1190 mm
Thickness folded	70 mm
Thickness opened	65 mm
Width	440 mm
Material	Aluminum
Loading capacity	150 kg
Weight	7,9 ± 0,4 kg

<sup>1</sup> Length measures are subjected to tolerances of ± 10 mm

Z 5000 3 BLADES PICK UP  
STRETCHER W/3 BELTS  
**ZS00400A**

Rev.2 (15/03/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>  
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222  
[www.spencer.it](http://www.spencer.it)