

Z2000

Foldable stretcher with body bag



Foldable stretcher for recovery of corpses, equipped with body bag.

Specific features

Resistant PVC fabric and easy to clean
 Footrest and wheels on the bottom, allow to move the device by tilting and sliding it on the ground
 The device can be folded, minimizing storage space
 System to lock the stretcher in open position easy to operate

Technical data

Length	2005 ± 10 mm
Width	500 ± 10 mm
Thickness w/footrest opened	420 ± 10 mm
Dimensions folded	1070x495x140±20 mm
Loading capacity	130 kg
Weight	9,5 ± 1 kg

Z2000 W/PVC GREY BODY BAG
 AND FOOTREST
ZS00500A

Rev.3 (14/06/2017)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

mortuary.spencer.it