

B-COVER



Cover for stretcher with basin.

Specific features

Lightweight and non-deformable
Four metal clips allow quick coupling with B-Support stretcher

Technical data

Length	1950 ± 10 mm
Width	570 ± 10 mm
Height	320 ± 10 mm
Material	HDPE

B-COVER – COVER FOR STRETCHER
WITH BASIN
ZS01000C

Rev.0 (02/11/2022)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222

mortuary.spencer.it