

Z 754

Head riser



System to support the head or limbs during autopsy or dressing procedures.

Specific Characteristics

High friction material to ensure great stability
 Easy to clean
 6 different radiuses for maximum versatility

Technical data

Dimensions	167x125x110 mm
Weight	730 g
Material	PU

Z 754 BACKREST LIFT
ZV00300B

Rev.1 (15/03/2018)

UNCHECKED COPY – further revisions will be available on <http://support.spencer.it>
 Spencer Italia S.r.l. Sala Baganza (PR) Italia Tel. +39.0521.541111 Fax +39.0521.541222
mortuary.spencer.it