

## Structural modification

### M SERIES

#### Internal pouches



CB00406A



CB00402A



CB00403A



CB00405A



CB00404A

The M Series internal pouches for bags and backpacks, also indicated for use as containers for other emergency material and first aid items, will no longer be equipped with transport handle/s.

The adjustment has been made in order to cancel the risk that the bags remain entangled in other items or rescue equipment during emergency operations and therefore reducing any consequent danger for patient and rescuer.

#### TECHNICAL DATA

Dimensions (M1, M2, M4, M5) (mm)	300x190xh100
Dimensions (M3) (mm)	300x190xh200
Weight (g) (M1, M2, M4, M5)	200
Weight (g) (M3)	550
Capacity (M1, M2, M4, M5) (L)	6
Capacity (M3) (L)	12
External material	Nylon Ripstop
Structural material	PE
Internal material Taffetà Nylon/ Cordura	Yes
Velcro® ridge fixation	Yes
Zipper closure	Yes
Zip pull	Yes
External transparent pocket	Yes
Compartments / internal pockets (M1)	12
Compartments / internal pockets (M2, M3)	7
Compartments / internal pockets (M4)	3
Compartments / internal pockets (M5)	6
Suitable backpack	Mark, R-aid
Elastic loops (M2)	5
Elastic loops (M3)	100
M 1	CB00406A
M 2	CB00402A
M 3	CB00403A
M 4	CB00405A
M 5	CB00404A

The revision described is valid from 27/02/2012