

Structural revision

MARK

Professional emergency rucksack



The revisions have been made on the inner parts of the rucksack . The product has undergone important structural and technical revision which offer important upgrades with respect to the previous style. Both the strength of the loading frame and comfort for the user have been improved and these aspects offer an important contribution to the functionality of the rucksack.

In specific the revisions are:

- Loading frame of the rucksack has been integrated in areas of contact
- Back part structured in lightweight aluminium alloy
- Padded shoulder straps
- Shoulder straps easily adjusted in height offering the possibility to adjust to meet ideal user comfort.
- Lower back area assisted with added padding
- High visibility reflective inserts
- Safety precaution labelling for safety of user and rucksack.

TECHNICAL DATAs

Dimensions	360 x 300 x h550 mm
Volume	59 L
Weight empty	3 kg
Loading capacity	Cordura
Structural material	PE / Aluminium
Handle material	Polyester
Base reinforcement	TPU
Outer pockets	3
Inner pockets	5 (transparent)
Fixing of gas bottle	Si (max. 2 l)
Inner pockets	Integrated CB00455A Accessory for CB00450A
Pouch fixing	Strap
CB00450A	Mark yellow / black without pouches
CB00455A	Mark yellow / black with pouches

The revisions will be available as from 01/03/2013