

# MINIZAINO

Compact BLS backpack  
Backpack replacement



Compact professional backpack equipped with devices most commonly used in BLS procedures.

<b>Changes</b>
<b>Backpack replacement</b>
The previous MINIZAINO is replaced by the A911 model
The included equipment is the same and is visible <a href="#">HERE</a>
Changes active for orders starting from 14/03/2016

<b>Specific Characteristics</b>
Dual transport possibility thanks to the padded shoulder straps and the soft carrying handle placed on the upper side
Adjustable shoulder straps with adjustable thoracic and abdominal cross bands with quick release buckle to provide more stability during transport
High visibility thanks to the reflective bands
The frontal pocket is capacious and can be completely opened. On both internal sides are sewn elastic bands to attach accessories
Large pocket outside the frontal flap equipped with zip closure
The main compartment is padded and is characterized by
- elastic bands on the perimeter suitable for accessories storage
- zip closure that allow the complete opening of the compartment
- strap bands to keep in place the internal pockets
- perimetral reinforcement to limit deformations
2 wide lateral pockets with zip closure
5 internal detachable pockets with zip closure and transparent top
Non-slip rubber coated bottom to make cleaning easier

<b>Technical data</b>	
Dimensions	490x470x210 ± 30 mm
Dimensions of each internal pocket	280x120x70 ± 10 mm
Materials	Nylon, TPU, PVC
Volume	34 ± 1 l
Weight of empty backpack w/pockets	1,97 ± 0,1 kg

Medical devices CE marked according to dir. 93/42/CE
A911 - COMPLETE BLS BACKPACK W/5 POCKETS
CB00913A